

## Auto Kerato-Refractometer





# The next generation of refractive care

# Function and aesthetics combined in the KR-1

The KR-1 is a revolutionary auto Kerato-Refractometer. The tiltable and rotatable touchscreen will change the way of working forever.

The operator can use the KR-1 from various angles, ensuring the best interaction with the patient. The Topcon KR-1 can be used in several set-ups, which use a minimum of working space.

#### **Features**

- Fully automated operation with touchscreen
- Flexible & space saving layout
- Reliable & fast measurements
- Compact ergonomic design







### Unique features of the KR-1




# KR-1 Fully automated with touchscreen

The KR-1 is fully operated through a 8.5 inch wide touchscreen. Simply touch the center of the pupil on the LCD screen to start measuring both eyes. The KR-1 automatically obtains objective refraction data of the right and left eye.




Fully automated! Just touch the center of the pupil




The KR-1 can detect automatically abnormal conditions of an eye, i.e. cataract. The instrument changes to 'cataract mode' and continues the measurement. The KR-1 incorporates a button for cataract, in case you would like to start measurements in the cataract mode from the start.

#### Manual mode

If required, the KR-1 can be set in manual mode as well. The manual mode is used for focusing only. The up/down and left/right movement will remain automated.





#### Easy to use color touchscreen control panel

All operating procedures can be performed with the touchscreen. The buttons on the screen are clear and easy to understand.



#### **KR-1** features



Wall set-up / behind patient

# 

Classic / conventional set-up

#### Flexible and space saving layout

The adjustable touchscreen control panel enables the operator to be positioned anywhere around the patient, because the control panel can be faced in a number of different directions. The KR-1 can be positioned in a conventional set-up, against a wall, or even in a corner of a room.

The compact body enables the operator to support the patient easily from many different positions during measurement even when the patient's eyelid needs to be opened. These unique aspects will contribute to space saving and flexible layout in your eye examination room.



Corner set-up



**Conventional position** 

(seated)

Front position (seated)




Ergonomic operation in standing position




Side position



Rear position





#### Reliable & fast measurement

With Topcon's exclusive Rotary Prism Technology, the KR-1 provides unmatched accuracy and reliability. This technology ensures reproducible results. The objective measurement of both eyes is done in a fast and comfortable way for the patient.


#### **Peripheral KRT**

The KR-1 allows the operator to take peripheral KRT measurements of the corneal curvature radius, in addition to inner KRT ring measurements, which is useful in corneal evaluation of contact lens fitting. The KR-1 can measure patients with a pupil as small as Ø2 mm.

#### Compact ergonomic design

The KR-1's new ergonomic design provides a more comfortable position for the patient with a 5° inclination of the patient's head, and optical head of the KR-1. The compact body of KR-1 and improved chinrest design enables easier access to patients. The operator will have better interaction with the patient.

Ergonomic optical head and headrest design



IMAGEnet



Compact body








#### Connectivity

The KR-1 can be connected directly to other instruments such as the Topcon CV-5000 automated phoropter. The KR-1 can be integrated into the Topcon Imaging and database software IMAGEnet. This software enables easy and fast access to patient data and is a useful tool to maintain customer contacts.

#### System chart





Subject to change in design and/or specifications without advanced notice. In order to obtain the best results with this instrument, please be sure to review all user instructions prior to operation.

Medical device class Im. Manufacturer: Topcon Corporation.







#### Topcon Europe Medical B.V.

Essebaan 11; 2908 LJ Capelle a/d IJssel; P.O. Box 145; 2900 AC Capelle a/d IJssel; The Netherlands Phone: +31-(0)10-4585077; Fax: +31-(0)10-4585045 E-mail: medical@topcon.eu; www.topcon-medical.eu

#### Topcon Danmark

Præstemarksvænge 25, 4000 Roskilde, Denmark Phone: +45-46-327500; Fax: +45-46-327555 E-mail: info@topcon.dk; www.topcon.dk

#### Topcon Scandinavia A.B.

Neongatan 2; P.O. Box 25; 43151 Mölndal, Sweden Phone: +46-(0)31-7109200; Fax: +46-(0)31-7109249 E-mail: medical@topcon.se; www.topcon.se

#### Topcon España S.A.

HEAD OFFICE; Frederic Mompou, 4; 08960 Sant Just Desvern; Barcelona, Spain Phone: +34-93-4734057; Fax: +34-93-4733932 E-mail: medica@topcon.es; www.topcon.es

#### Topcon Italy

Viale dell' Industria 60; 20037 Paderno Dugnano, (MI) Italy Phone: +39-02-9186671; Fax: +39-02-91081091 E-mail: info@topcon.it; www.topcon.it

#### Topcon France Medical S.A.S.

BAT AI; 3 Route de la Révolte, 93206 Saint Denis Cedex Phone: +33-(0)1-49212323; Fax: +33-(0)1-49212324 E-mail: topcon@topcon.fr; www.topcon-medical.fr

#### Topcon Deutschland GmbH

Hanns-Martin-Schleyer Strasse 41; D-47877 Willich, Germany Phone: (+49) 2154-885-0; Fax: (+49) 2154-885-177 E-mail: info@topcon-medical.de; www.topcon-medical.de

#### Topcon Polska Sp. z o.o.

ul. Warszawska 23; 42-470 Siewierz; Poland Phone: +48-(0)32-670-50-45; Fax: +48-(0)32-671-34-05 www.topcon-polska.pl

#### **Topcon Great Britain Medical Ltd.**Topcon House; Kennet Side; Bone Lane; Newbury

Phone: +44-(0)1635-551120; Fax: +44-(0)1635-551170

E-mail: medical@topcon.co.uk, www.topcon.co.uk

#### **Topcon Ireland**

Unit 276, Blanchardstown; Corporate Park 2 Ballycoolin; Dublin 15, Ireland Phone: +353-18975900; Fax: +353-18293915 E-mail: medical@topcon.ie; www.topcon.ie

<sup>\* -25</sup>D < spherical refractive power + cylindrical refractive power or spherical refractive power + cylindrical refractive power < +22D