Outstanding Image Quality in a Fraction of the Time

Leading retina specialists share their impressions of PixelSmart[™], the latest image processing algorithm from Topcon Healthcare

Sponsored content by Topcon

Topcon Healthcare, an Optical Coherence Tomography (OCT) instrumentation leader, offers both spectral domain (SD) and sweptsource (SS) OCT options, has introduced a new image processing algorithm to its unique DRI OCT Triton[™] series. The novel PixelSmart technology takes SS-OCT imaging to the next level and offers unparalleled image quality by reducing speckle noise and improving image contrast. Many OCT devices repeatedly scan each retinal location and

average the data, in order to obtain

satisfactory image quality. Clinicians using such devices will be excited to try PixelSmart which can provide a detailed assessment of the posterior segment both quickly and efficiently. Thanks to Triton's exceptional highdensity SS-OCT data, the new image processing algorithm is capable of producing rich, detailed images without compromising speed or the scan area.

Here, ophthalmology experts specializing in retinal care present clinical cases assessed with the Triton utilizing PixelSmart.

Case One

High myopia with choroidal neovascular membrane

Luis Arias, MD, PhD, Head of the Retina Department of University Hospital of Bellvitge and Aggregate Professor of Ophthalmology, University of Barcelona, Spain, presents a case of a 66-year-old male patient suffering from high myopia and a choroidal neovascular membrane, with a previous history of photodynamic treatment, 20 years prior.

Professor Arias shares his view on the utility of the new image processing technology available from Topcon Healthcare, "With a good compromise between acquisition time and image quality, PixelSmart can improve workflow in busy clinics. Within a couple of seconds, I can analyze high density scans with good image quality, allowing me to evaluate details of the vitreous, retina and choroid."



Figure I. High myopia and choroidal neovascular membrane. Triton Fundus image captured simultaneously with the OCT.









Figure 2. High myopia and choroidal neovascular membrane A) Triton 3D OCT scan, raw image, B) Triton 3D OCT scan, PixelSmart image, C) Triton OCT radial scan, averaged image.

topconhealthcare.com

Case Two

Diabetic macular edema

José Maria Ruiz-Moreno, MD, PhD, Head of Ophthalmology, Puerta de Hierro University Hospital, Majadahonda, Madrid, Professor at the University of Castilla-La Mancha Medical School, Albacete, Spain, presents a case of a 77-year-old female patient with diabetic macular edema.

On the novel imaging algorithm, Professor Moreno comments, "PixelSmart is a great tool for screening the retina. It improves the image quality of standard 3D scans, aiding in the assessment of retinal diseases."













Figure 3. Diabetic macular edema F) Triton 3D volume scan, raw image, G) Triton 3D volume scan, PixelSmart image. Averaged line scans do not allow 3D reconstruction - it is an advantage of 3D cube scans.

Figure 3. Diabetic macular edema A) Triton 3D OCT scan, raw image, B) Triton 3D OCT scan, PixelSmart image, C) Triton OCT line scan, averaged image.







Figure 3. Diabetic macular edema D) Triton 3D scan, Enface raw image, E) Triton 3D scan, Enface PixelSmart image. Deep vascular complex segmentation (IPL/INL 15.6-IPL/INL 70.2).

topconhealthcare.com

Case Three

Wet age-related macular degeneration

Heloísa Nascimento, MD, PhD, from the Department of Ophthalmology, Federal University of São Paulo, Paulista Medical School – UNIFESP/EPM, Brazil, presents a case study of an 80-yearold female patient with wet age-related macular degeneration.

Figure 4. Wet age-related macular degeneration. Triton Fundus image captured simultaneously with the OCT.







Figure 5. Wet age-related macular degeneration A) Triton 3D OCT scan, raw image, B) Triton 3D OCT scan, PixelSmart image.

Heloísa Nascimento, MD, PhD, also shares a case study of a 33-year-old male patient diagnosed with ocular syphilis. Professor Nascimento shares her

impressions of Topcon Heathcare's new image processing technology, "PixelSmart has greatly improved the image quality of 3D scans. Having good images with a high-density scan allows for a better evaluation of the retinal and choroidal alterations. This is especially important to image patients who do not collaborate with the exam or have poor fixation, when it is difficult to position a line or radial high-quality scan exactly where the alteration is."

Blink-of-an-eye imaging

Case Four

Ocular syphilis

With PixelSmart technology, Topcon Healthcare once again pushes the boundaries of OCT imaging, providing unprecedented image quality available to clinicians for instant, precise analysis, improving the clinic's performance and ultimately – patients' outcomes. It further cements Topcon Healthcare's position as the market leader in the field of OCT, constantly innovating and improving its portfolio to deliver the best solutions to ophthalmic practices around the world.

Topcon Healthcare recommends installing the viewing system with PixelSmart on a network to allow image review from exams rooms in addition to the imaging room.

Not all products, services or offers are approved or offered in every market and products vary from one country to another. Contact your local distributor for countryspecific information.







Figure 6. Ocular syphilis. Triton Fundus image captured simultaneously with the OCT.



Figure 7. Ocular syphilis A) Triton 3D OCT scan, raw image, B) Triton 3D OCT scan, PixelSmart image.

topconhealthcare.com