## DRI OCT Triton™ Series

Swept Source Optical Coherence Tomography



**NEW** PixelSmart™



## See. Discover. Explore.

The diagnostic power of Swept Source OCT Deep Range Imaging<sup>1</sup>.



"Swept Source adds a new dimension to OCT. The TOPCON DRI Swept Source OCT is easy to use, provides unique clinical information, and has improved my practice. For the first time, we can in-vivo visualize not only the vitreo-retinal interface but also the cortical vitreous which is important at the time when more and more therapies are delivered via intra-vitreal injections. Deeper imaging brings choroidal thickness, helping guide my clinical decisions. Seeing more helps guide my therapy and allows me to treat more effectively. I find Swept Source OCT an essential tool to look for biomarkers of disease regression or progression."

#### Prof. P. E. Stanga

Manchester Royal Eye Hospital, Manchester Vision Regeneration (MVR) Lab at NIHR/ Welcome Trust Manchester CRF & University of Manchester

1) Fabio Lavinsky, Daniel Lavinsky. Novel perspectives on swept-source optical coherence tomography. Int J Retin Vitr (2016) 2:25

## Welcome to the New Frontier in OCT Imaging.

The DRI OCT Triton combines the world's first<sup>2</sup> Swept Source OCT technology with multimodal fundus imaging.

#### **Image Quality**

Triton's Swept Source with its 100 kHz scanning speed and 1,050nm wavelength results in a clear and detailed images even for the deepest layers of the eye with short acquisition time. Visualize not only the retina and vitreous, but also the choroid and sclera<sup>1</sup>.

#### Diagnostic capability

Seeing deeper makes it possible to have a better understanding of many ocular pathologies<sup>1</sup>. With features such as OCT angiography, fundus autofluorescence and en face OCT, Triton empowers clinicians with multimodal imaging capability to help assess and preserve patient's eye health.


#### **Practice efficiency**

The Triton's automated functions, including single scan captures and SMARTTrack™ system, are designed to optimize your practice workflow by simplifying data capture, analysis and follow-up.




## See Deeper<sup>1</sup>. See More.

#### **Proliferative diabetic retinopathy**




Courtesy: Prof. P. E. Stanga, Manchester Royal Eye Hospital, Manchester Vision Regeneration (MVR) Lab at N IHR/Welcome Trust Manchester CRF & University of Manchester




Courtesy: Prof. P. E. Stanga, Manchester Royal Eye Hospital, Manchester Vision Regeneration (MVR) Lab at N IHR/Welcome Trust Manchester CRF & University of Manchester

#### **Central serous chorioretinopathy**



Courtesy: Prof. P. E. Stanga, Manchester Royal Eye Hospital, Manchester Vision Regeneration (MVR) Lab at N IHR/Welcome Trust Manchester CRF & University of Manchester



Courtesy: Prof. P. E. Stanga, Manchester Royal Eye Hospital, Manchester Vision Regeneration (MVR) Lab at N IHR/Welcome Trust Manchester CRF & University of Manchester


4

<sup>\*</sup> FA photography and FAF photography can be performed using only DRI OCT Triton Plus.

<sup>\*</sup> FA photography and FAF photography can be performed using only DRI OCT Triton Plus.


## See Deeper¹. See More.

#### Pathological myopia




Courtesy: Prof. P. E. Stanga, Manchester Royal Eye Hospital, Manchester Vision Regeneration (MVR) Lab at N IHR/Welcome Trust Manchester CRF & University of Manchester

#### Macular pucker




Courtesy: Prof. P. E. Stanga, Manchester Royal Eye Hospital, Manchester Vision Regeneration (MVR) Lab at N IHR/Welcome Trust Manchester CRF & University of Manchester

#### Image through cataract



a)



a, b, c courtesy of Kazuya Yamagishi, MD (Hirakata Yamagishi Eye Clinic, Japan)

1) Fabio Lavinsky, Daniel Lavinsky, Novel perspectives on swept-source optical coherence tomography. Int J Retin Vitr (2016) 2:25

## Discover What Lies Beneath

TOPCON's SS OCT Angio™ combines OCT angiography with a Swept Source OCT.

OCTARA™, a proprietary image processing algorithm, provides highly sensitive angiographic detection³, allowing for visualization of vascular structures even in the choroid and deeper retinal layers.

#### High-sensitivity Imaging and Deeper Intravascular Flow Visualization<sup>1</sup>

Swept Source technology and OCTARA<sup>TM</sup> allow the deeper structures to be visualized with less depth-dependent signal roll-off<sup>3</sup>. Additionally, the 1 $\mu$ m wavelength makes OCT Angiography imaging possible for patients with media opacities.

#### Rapid Scanning, Real Time Eye Tracking

At 100,000 A-Scans per second coupled with invisible\* scanning lines and the SMARTTrack™ eye tracking system, 'the Triton quickly captures a dense data set and provides an en face OCT Angiography image of the retinal microvascular flow network'<sup>3</sup>.


\* OCT Angiography scanning line may be faintly visible during capture to some people with certain conditions

#### **Efficiency & Workflow Integration**


Multimodal platform provides comparison of microvascular impairment with FA, FAF, OCT and true color fundus images in a single device\*.

\*DRI OCT Triton Plus

#### 12x12mm 512 pixels




Courtesy: Akiniro Ishibazawa, MD, PhD, Asahikawa Medical Universil Graduate School of Medical Sciences, Hokkaido, Japan





o, Japan Retinal Research and Imaging, Austria

#### OCTARA™

OCTARA<sup>TM</sup> is the image processing technology which extracts the signal changes derived from vascular flow using multiple OCT B-Scans acquired at the same position. It demonstrates high sensitivity for the detection of low blood flow in microvasculature $^3$ .








ourtesy: SriniVas R. Sadda, M.D., Doheny Eye Institute, UCLA

Courtesy: SriniVas R. Sadda, M.D., Doheny Eye Institute, UCLA

Courtesy: Kazuya Yamagishi, MD., Hirakata Yamagishi Eye Clinic, Jap

#### **Multimodal Viewing**

En face angiography images, B-scans and fundus photography can all be viewed on a single screen using IMAGEnet®6 and PinPoint™ registration, so that area of interest can be assessed using multiple image modalities. Selected layers can easily be customized to enhance the clarity of specific pathological features.




<sup>3)</sup> Magdy Moussa, Mahmoud Leila, Hagar Khalid. Imaging choroidal neovascular membrane using en face swept-source optical coherence tomography angiography. Clinical Ophthalmology 2017:11 1859-1869

<sup>3)</sup> Magdy Moussa, Mahmoud Leila, Hagar Khalid. Imaging choroidal neovascular membrane using en face swept-source optical coherence tomography angiography. Clinical Ophthalmology 2017:11 1859–1869


## Discover What Lies Beneath

#### **Proliferative diabetic retinopathy**




SS OCT Angio™ Montage
Courtesy: Akihiro Ishibazawa, MD, PhD.
Asahikawa Medical University Graduate School of Medical Sciences, Hokkaido, Japan

#### Before treatment




Courtesy: Akihiro Ishibazawa, MD, PhD. Asahikawa Medical University Graduate School of Medical Sciences, Hokkaido, Japan

#### After treatment



Courtesy: Akihiro Ishibazawa, MD, PhD. Asahikawa Medical University Graduate School of Medical Sciences, Hokkaido, Japan

#### **Branch retinal vein occlusion**





Courtesy:Yuichiro Ogura, MD, Professor and Chairman, Department of Ophthalmology and Visual Science, Nagoya City University, Nagoya, Japan


10 11

## Discover more possibilities: see beyond and deeper

OCTA metrics on Triton SS OCT Angio allows clinicians to objectively and quantitatively assess retinal vasculature, providing valuable insights into the patient's eye health.










## Swept Source OCT Imaging

#### 1,050nm wavelength

The longer wavelength light provides better tissue penetration, allowing visualization into the deepest layers of the eye1.






Courtesy: Professor Jose Maria Ruiz Moreno MD, University of Albacete, Spain


#### **Swept Source OCT technology;** scanning speed of 100,000 A-scans/sec

The fast scanning speed of 100,000 A-scans/sec enables capture of clear B-scans<sup>4</sup> by acquiring more A-scans within a given image acquisition time. This helps to reduce artifacts from involuntary eye movements such as saccades and blinks.



#### Invisible scan lines

The invisible 1,050nm wavelength light helps patients concentrate on the fixation target during the scan, reducing involuntary eye movement. It supports more efficient workflow in a practice by reducing the need to rescan.




<sup>1)</sup> Fabio Lavinsky, Daniel Lavinsky, Novel perspectives on swept-source optical coherence tomography. Int J Retin Vitr (2016) 2:25

<sup>4)</sup> Shoji Kishi. Impact of swept source optical coherence tomography on ophthalmology. Taiwan Journal of Ophthalmology 6 (2016) 58-68


#### **En face OCT imaging**

en face imaging allows for independent dissection of a depth range defined by two boundaries, selected from seven possible boundaries, by flattening the 3D data cube.





Courtesy: Prof. T. Nakazawa, MD,PhD, Tohoku University, Japan




Courtesy: Prof. T. Nakazawa, MD,PhD, Tohoku University, Japan

#### Vitreous visualization

#### Dynamic Focus™

Dynamic Focus<sup>™</sup> on Triton allows for acquisition of images with near uniform focus and image quality throughout the entire depth of the image, which helps to enhance the typically weaker signal from the vitreous.







#### **EVV (Enhanced Vitreous Visualization™)**

EVV helps clinicians assess vitreous and vitreoretinal interface abnormalities<sup>1</sup>. Contrast can be quickly adjusted to the needs of the physician, depending on the area of greatest interest.





## Discover from Cornea to Choroid



Courtesy: Michael H. Chen, O.D.

#### **Anterior segment imaging**


Triton's (optional) anterior segment imaging capabilities allow for visualization of the cornea, anterior chamber angle, iris and sclera<sup>5</sup>.

#### Image samples


#### OCT image B-scan length 16mm



#### Anterior segment in Radial scan




#### Anterior segment in 3D scan



#### 7 boundary segmentation/5 layer thickness maps/caliper function

Retinal tissue layers are automatically segmented by the Topcon Advanced Boundary Software (TABS™), enabling quantification of retinal thickness and sub layers<sup>67</sup>.




#### **Choroidal thickness maps**

Understanding choroidal structural changes in relation to various disease and disease management, has been of long-standing interest to clinicians<sup>8</sup>. Triton provides clear visualization of the choroid, and generates choroidal thickness maps to aid in the assessment of choroidal structure.

Retina	between the ILM-OS/RPE boundaries
RNFL	between the ILM-RNFL/GCL boundaries
GCL+	between the RNFL/GCL-IPL/INL boundaries
GCL++	between the ILM-IPL/INL boundaries
CSI	between the BM-CSI boundaries or ILM-CSI boundaries

#### **Hood Report (for Glaucoma)**


Retinal Thickness/RNFL/GCL and Optic Nerve Metrics in just one scan. This report streamlines the decision-making process through the correlation of structure (GCC/RNFL) with function (overlay of visual field test locations)<sup>6</sup>.



## DRI Meets Multimodal Fundus Imaging: See the Whole Picture


#### Swept Source OCT incorporates multimodal fundus imaging

DRI OCT Triton can acquire the OCT and fundus image in a single capture.
PinPoint™ Registration identifies the location of the B-scan on the fundus image. Comparison between the B-scan and fundus image can support clinical efficiency during diagnosis.



Courtesy: Jay M. Haynie, O.D.

OCT + Color fundus



Courtesy: Jay M. Haynie, O.D.

OCT + FAF

#### True color\* Fundus images

The DRI OCT Triton offers a true color, non-mydriatic fundus image. Fluorescein Angiography (FA) and Fundus Autofluorescence (FAF) are available\*\* to enhance the diagnostic capability of Triton Plus. The all-in-one device supports efficient workflow in practice.


\*Color fundus image with white light, with 24-bit color.
\*\*DRI OCT Triton Plus:

OCT /Anterior OCT (Option)/ OCT Angiography (Option) /Color / Red-Free / FA / FAF DRI OCT Triton :

OCT /Anterior OCT (Option)/ OCT Angiography (Option) /Color/ Red-Free

#### **Sensor Information**

Color Imege Sensor: 5MP FA Image Sensor: 4MP FAF Image Sensor: 4MP



Zhichao Wu, Denis S. D. Weng, Rashmi Rajshekhar, Abinaya Thenappan, Robert Ritch, Donald C. Hood. Evaluation of a Qualitative Approach for Detecting Glaucomatous Progression Using Wide-Field Optical Coherence Tomography Scans. Trans Vis Sci Tech. 2018;7(3):5.
 Beatriz Abadia, Ines Suñen, Pilar Calvo, Francisco Bartol, Guayente Verdes, Antonio Ferreras. Choroidal thickness measured using swept-source optical coherence

<sup>7)</sup> Beatriz Abadia, Ines Suñen, Pilar Calvo, Francisco Bartol, Guayente Verdes, Antonio Ferreras. Choroidal thickness measured using swept-source optical coherence tomography is reduced in patients with type 2 diabetes. PLoS ONE 13(2): e0191977.

<sup>8)</sup> Sushmitha Rao Uppugunduri, Mohammed Abdul Rasheed, Ashutosh Richhariya et al. Automated quantification of Haller's layer in choroid using swept-source optical coherence tomography. PLoS ONE 13(3):e0193324

## Smarter Tracking. Smarter Workflow.

#### $\mathbf{SMARTTrack}^{\!\scriptscriptstyle{\mathsf{M}}}$

The SMARTTrack™ system enhances the tracking and follow-up ability of the Triton with a variety of functions designed to enhance its user-friendliness:

- Fundus-Guided Acquisition (FGA)
- Follow-up Function
- Tracking Photography


#### Stereo photography

Three dimensional visualization of color fundus images can be achieved by acquiring images in stereo photography mode. Triton's on-screen acquisition guidance supports quick and easy operation with auto alignment for capturing stereo pairs.





#### Panoramic wide field photography

In addition to macular and disc imaging, the Triton provides wide coverage of the retina. A panoramic graphic can be created from multiple fundus or OCT Angiography images.











#### **Fundus Guided Acquisition (FGA)**

OCT scan location can be easily set by selecting the scan area on the fundus image. With FGA, the operator can choose to capture or import a fundus image, select the scan location and automatically acquire a B-scan at that location.



#### **Follow-up Function**

This function allows you to retrieve and re-analyze the same location at follow-up, for comparison of past and current images. All an operator needs to do is simply select the past data and Triton automatically captures the same area.





#### Motion Correction/Compensation/ Rescanning Function

#### **Motion Correction**

Corrects the Z direction movement

#### **Compensation Function**

Tracks the eye and then compensates for the X direction movement.

#### **Rescanning Function**

Parts of the scanning area may be missed due to Y direction eye movement. In such a case, the rescanning function automatically activates to rescan the missing area.

# → **→**

Before compensation


After compensation


#### Alignment navigation

When an operator wishes to acquire an image, the Triton's monitor guides the operator to reduce potential errors and make operation simple.

- Auto focus and auto shoot, in color/FAF mode
- Auto focus, auto-Z and Z-lock function, in OCT mode







close

OK

too far

#### **Small pupil solution**

#### **Live Fundus View**

The fast scanning speed allows the Triton to create a live en face fundus image, an ideal tool for precisely visualizing the scan position. This makes the disc, retinal vessels and scanning position easy to see, even in patients with small pupils.

#### OCT capture mode without retinal photography

Triton can also capture a 3D scan, with or without color fundus photography, to avoid a miotic response and better meet the needs of patients with the smallest pupils.



20

# Reinventing Swept Source Imaging:

Triton with PixelSmart<sup>™</sup> is the next stage in Swept Source OCT imaging


#### **NEW! PixelSmart™**


Triton with PixelSmart
The next level in Swept Source OCT Imaging

PixelSmart is Topcon's new image processing algorithm which reduces speckle noise to improve contrast, producing a better image quality.

PixelSmart is available for all Triton 3D scans, existing and new:

- 3D Wide
- 3D Macula
- 3D Disc
- Combination scans





## Specifications

Observation & Photography of Fundus	
Photography Type	Color, FA*, FAF*, Red-free**
Picture Angle	45°
	Equivalent 30° (Digital Zoom )
Operating Distance	34.8mm
Photographable Diameter of Pupil	Normal: φ4.0mm or more
	Small pupil mode diameter: φ3.3mm or more
Observation & Photography of Fundus	Tomogram
Scanning Range (on fundus)	Horizontal Within 3 to 12mm
	Vertical Within 3 to 12mm
Scan Pattern	3D scan
	Linear scan (Line-scan/Cross-scan/Radial-scan/Raster-scan)
Scan Speed	100,000 A-Scans per second
Lateral Resolution	20μm
In-depth Resolution	Digital: 2.6µm
	Optical function: 8µm
Photographable Diameter of Pupil	φ2.5mm or more
Observation & Photography of Fundus	Image / Fundus Tomogram
Fixation target	Internal fixation target :
	Dot matrix type organic EL
	The display position can be changed and adjusted.
	The displaying format can be changed.
	Peripheral fixation target :
	This is displayed according to the internal fixation target
	indicated position.
	External fixation target
Observation & photography of anterior	-
Photography type	IR
Operating distance	17mm
Observation & photography of anterior	11.11.11.11
Operating distance	17mm
Scan range (on cornea)	Horizontal Within 3 to 16mm
Scarrange (on comea)	Vertical Within 3 to 16mm
Scan pattern	3D scan
Scall pattern	
Coop speed	Linear scan (Line-scan/Radial-scan)
Scan speed	100,000 A-Scans per second
Fixation target	Internal fixation target
Electric Debies	External fixation target
Electric Rating	14 ly 100 0 10 V
Power Source	Voltage: 100-240V
	Frequency: 50-60Hz
Power input	250VA
Dimensions / Weight	
Dimensions	320-359 mm(W) X 523-554 mm(D) X 560-590 mm(H)
Weight	21.8kg (DRI OCT Triton)
	23.8kg(DRI OCT Triton Plus)

 $<sup>^{\</sup>ast}$   $\,$  FA photography and FAF photography can be performed only with the DRI OCT Triton Plus.

<sup>\*\*</sup> In this digital red-free photography, the color image is processed and is displayed as a pseudo red-free photographed image.

\*\*\* Observation & photography of anterior segment can be performed only when the anterior segment attachment is used.











#### TOPCON CORPORATION

75-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8580, JAPAN. Phone: +81-(0)3-3558-2522/2502 Fax: +81-(0)3-3965-6898 www.topconhealthcare.jp

#### TOPCON HEALTHCARE SOLUTIONS EMEA OY TOPCON SCANDINAVIA

Saaristonkatu 9, 90100 Oulu, FINLAND Phone: +358-20-734-8190 E-mail: thsemea.sales@topcon.com www.topconhealthcare.eu

### © RD TOPCON EUROPE MEDICAL B.V. Essebaan 11, 2908 LJ Capelle a/d IJssel P.O. Box 145, 2900 AC Capelle a/d IJssel THE NETHERLANDS

Phone: +31 -(0)10-4585077 Fax: +31 -(0)10-4585045 E-mail: medical@topcon.com www.topconhealthcare.eu

#### TOPCON DANMARK

OPCON DANMARK
Praestemarksvej 25, 4000 Roskilde
DANMARK
Phone: +45-46-327500 Fax: +45-46-327555 E-mail: info.todk@topcon.com www.topconhealthcare.eu

Neongatan 2, P.O.Box 25, 43151 Mölndal SWEDEN Phone: +46-(0)31-7109200 Fax: +46-(0)31-7109249 E-mail: info.hcs.se@topcon.com www.topconhealthcare.eu

#### TOPCON ESPAÑA

Oficina Principal en España Frederic Mompou, 4, 08960 Sant Just Desvern, Barcelona, SPAIN Phone: +34-93-4734057 Fax: +34-93-4733932 E-mail: medica@topcon.com www.topconhealthcare.eu/es\_ES

#### TOPCON ITALY

Viale dell' Industria 60. 20037 Paderno Dugnano, (MI) ITALY Phone: +39-02-9186671 Fax: +39-02-91081091 E-mail: info@topcon.com www.topconhealthcare.eu

#### TOPCON FRANCE MEDICAL

1 rue des Vergers, Parc Swen, Bâtiment 2, 69760 Limonest, FRANCE Phone: +33-(0)4-37 58 19 40 Fax: +33-(0)4-72 23 86 60 E-mail: topconfrance@topcon.com www.topconhealthcare.eu

#### TOPCON DEUTSCHLAND MEDICAL

Hanns-Martin-Schleyer Strasse 41, D-47877 Willich, GERMANY Phone: (+49)2154-885-0 Fax: (+49)2154-885-177 E-mail: info@topcon-medical.de www.topconhealthcare.eu

#### TOPCON POLSKA SP. Z. O. O.

ul. Warszawska 23, 42-470 Siewierz POLAND Phone: +48-(0)32-670-50-45 Fax: +48-(0)32-671-34-05 E-mail: info.tpl@topcon.com www.topconhealthcare.eu

#### TOPCON (GREAT BRITAIN) MEDICAL

Topcon House, Kennet Side, Bone Lane, Newbury, Berkshire RG14 5PX UNITED KINGDOM Phone: +44-(0)1635-551120
Fax: +44-(0)1635-551170
E-mail: medical.tgbm@topcon.com
www.topconhealthcare.eu

#### TOPCON IRELAND MEDICAL

Unit 292, Block G, Blanchardstown, Corporate Park 2 Ballycoolin Dublin 15, D15 DX58, IRELAND Phone: +353-12233280 E-mail: medical.ie@topcon.com www.topconhealthcare.eu

IMPORTANT In order to obtain the best results with this instrument, please be sure to review all user instructions prior to operation.

