

Henson 9000

Technical Specifications



Henson 9000 Perimeter: Technical specifications

Test specifications

Visual field test range	60° (monocular) / 160° (binocular)
Visual field testing distance	25 cm
Stimulus intensity (maximum)	10,000 asb
Background illumination	31.4 asb
Stimulus duration	200 ms
Stimulus size	Goldmann III
Stimulus colour	White
Test method	Standard Automated Perimetry (SAP)

Screening tests/patterns

Smart Supra - single stimulus	30, 64 and 86 point tests
Smart Supra - multiple stimulus	30, 64 and 86 point tests
Esterman (Driving)	Groups 1 (120 point) and 2 (124 point) (EU standard)
Customised tests	Test locations can be manually added to all Smart Supra screening tests

Threshold tests/patterns

ZATA Standard - threshold central	10-48**; 24/30-2 (extendable in-test)
ZATA Fast - threshold central	10-48; 24/30-2 (extendable in-test)

Average test times*

Smart Supra - single stimulus	-1 min (30 points); -3.5 minutes for fully extended 86 point test
Smart Supra - multiple stimulus	Under 30 seconds (30 points)
ZATA (24-2)	- 3 minutes per eye
ZATA Fast (24-2)	- 2.5 minutes per eye

Fixation control

Fixation target	Single or 4-point LED diamond pattern
Heijl-Krakau	Yes (ZATA)
Video eye monitor	Yes

Software features

Patient management database	MS Windows [®] compatible; networkable
Practice management integration	EMR compatibility (parameter passing and text file)
Hemifield analysis	Yes
Progression analysis	Yes
HFA data import	Yes

Connectivity

DICOM	Yes (images)
Ethernet	Yes, via connected computer
Database backup	Removable, network or cloud storage

Dimensions

Weight	13.5 kg
Size, W x D x H	440 x 400 x 452 (mm)

Classification

Medical device	Class I
Applied part	Type B
Control device	External PC / laptop / tablet running MS Windows [®] Professional, v.8, and above
Patient unit inputs/outputs	CI3 mains input; Patient Response Button; 2 x USB Type B connector
Electrical requirements	85 - 263 V AC, 50/60 Hz, 60 VA
Optional printer	Any compatible with controlling computer

* Test times are approximations only and will differ according to patient response time and level of loss.

**The 10-48 pattern presents 48 stimuli on a 3 degree square matrix, displaced 3 degrees from the horizontal and vertical mid-lines, within the central 10 degrees.

©Topcon Healthcare. All rights reserved.

Microsoft[®] and Windows[®] are either registered trademarks or trademarks of Microsoft Corporation[®] in the United States and/or other countries.

Not all products, services or offers are approved or offered in every market, and products vary from one country to another.

Contact your local distributor for country-specific information and availability.

TOPCON CORPORATION

75-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8580, JAPAN.
Phone: +81-(0)3-3558-2522/2502 Fax: +81-(0)3-3965-6898
www.topconhealthcare.jp

TOPCON HEALTHCARE SOLUTIONS EMEA OY

Saaristonkatu 9, 90100 Oulu, FINLAND
Phone: +358-20-734-8190
E-mail: thsemea.sales@topcon.com
www.topconhealthcare.eu

TOPCON EUROPE MEDICAL B.V.

Essebaan 11, 2908 LJ Capelle a/d IJssel
P.O. Box 145, 2900 AC Capelle a/d IJssel
THE NETHERLANDS
Phone: +31 -(0)10-4585077
Fax: +31 -(0)10-4585045
E-mail: medical@topcon.com
www.topconhealthcare.eu

TOPCON DANMARK

Praestemarksvej 25, 4000 Roskilde
DANMARK
Phone: +45-46-327500
Fax: +45-46-327555
E-mail: info.todk@topcon.com
www.topconhealthcare.eu

TOPCON SCANDINAVIA

Neongatan 2, P.O.Box 25,
43151 Mölndal, SWEDEN
Phone: +46-(0)31-7109200
Fax: +46-(0)31-7109249
E-mail: info.hcs.se@topcon.com
www.topconhealthcare.eu

TOPCON ESPAÑA

Oficina Principal en España
Frederic Mompou, 4,
08960 Sant Just Desvern, Barcelona, SPAIN
Phone: +34-93-4734057
Fax: +34-93-4733932
E-mail: medica@topcon.com
http://www.topconhealthcare.eu/es_ES

TOPCON ITALY

Viale dell' Industria 60,
20037 Paderno Dugnano, (MI) ITALY
Phone: +39-02-9186671
Fax: +39-02-91081091
E-mail: info@topcon.com
www.topconhealthcare.eu

VISIA IMAGING S.R.L.

Via Martiri della Libertà 95/e, 52027 San Giovanni Valdarno (AR), ITALY

TOPCON FRANCE MEDICAL

1 rue des Vergers, Parc Swen,
Bâtiment 2, 69760 Limonest, FRANCE
Phone: +33-(0)4-37 58 19 40
Fax: +33-(0)4-72 23 86 60
E-mail: topconfrance@topcon.com
www.topconhealthcare.eu

TOPCON DEUTSCHLAND MEDICAL

Hanns-Martin-Schleyer Strasse 41,
D-47877 Willich, GERMANY
Phone: (+49)2154-885-0
Fax: (+49)2154-885-177
E-mail: info@topcon-medical.de
www.topconhealthcare.eu

TOPCON POLSKA SP. Z. O. O.

ul. Warszawska 23, 42-470 Siewierz
POLAND
Phone: +48-(0)32-670-50-45
Fax: +48-(0)32-671-34-05
E-mail: info.tpl@topcon.com
www.topconhealthcare.eu

TOPCON (GREAT BRITAIN) MEDICAL

Topcon House, Kennet Side, Bone Lane,
Newbury, Berkshire RG14 5PX
UNITED KINGDOM
Phone: +44-(0)1635-551120
Fax: +44-(0)1635-551170
E-mail: medical.tgbm@topcon.com
www.topconhealthcare.eu

TOPCON IRELAND MEDICAL

Unit 292, Block G, Blanchardstown,
Corporate Park 2 Ballycoolin
Dublin 15, D15 DX58, IRELAND
Phone: +353-12233280
E-mail: medical.ie@topcon.com
www.topconhealthcare.eu